

ABSTRACT

Plasma-density interferometers are typically -

- robustly mounted, two-arm, high-maintenance devices

A **Second-Harmonic, Dispersion Interferometer** measures Df using 2w common-path beams



The 1-D SHDI uses CW (Nd:YAG, or CO₂) lasers ¹
A pulsed, high-intensity laser provides for large-area 2-D time-resolved imaging ^{2,3}

- 0.01-radian sensitivity, 0.1-mm resolution, 1-ns sampling, 10⁵ Hz


The 2-D SHDI is being characterized on pulsed-gas jets/plasma jets

- 2 x 3 cm, 50-us, n·dl ~ 10¹⁶ cm⁻²

DATA AQISITION



PLASMA SOURCE




Parker Series 9 Valve - 2A, 500 us
N_{gas} ~ 10¹⁷ cm⁻³
3-D printed insulator
Brass/SS electrodes
1-3 kV, 50 us, 500J square pulse
N_{plasma} ~ 10¹⁶ cm⁻³


Plasma-Density Imaging – 2-Dimensional, - Second-Harmonic, Dispersion Interferometer

C. T. Chavez, A. Egly, F. J. Wessel
L-Egant Solutions, LLC, Irvine, California, USA

DATA CAPTURE




Δφ Sample - Δφ Background
8.14mm




$$\Delta\phi(x, y) = 2.82 \times 10^{-13} \{ \lambda_{532} - \lambda_{1064} \} \int_{-z_0}^{z_0} N_e(x, y, z) dz.$$


15.9mm



PyAbel TRANSFORMS




GAS JET IN AIR



Gladstone-Dale Eqn.:

$$\bar{N} = \frac{\lambda}{4\pi L \Delta n_0} \Delta\phi$$

N(r, z)
Projection



Density Profile
N(x, y)
Peak Density:
1.18E19/cm³

CONCLUSIONS

- Scalable for measurements on plasma, gas, fluids
- Simple design, with a single laser source and minimal optics
- Common path 2w beams eliminates noise
- Integrated S/W processing facilitates real-time data analysis
- Results are much easier to interpret compared to standard “fringe counting”
- Simple design
- Robustly stable
- Low-maintenance
- Low-cost
- US Patent #11,221,293

REFERENCES

- 1 Brandi, F., Wessel, F.J., Lohff, C., Duff, J.R., Haralson, Z.O., 2020, App.Optics 59, p.8486-8493.
- 2Brandi, F. and Wessel, F.J., 2020, Optics Letters 45, p.4304.
- 3Wessel, F. and Brandi, F. US Patent #11221293.

ACKNOWLEDGEMENTS

Cameron Chavez: Cameron@l-egantsolutions.com
Frank Wessel: Frank@l-egantsolutions.com
Supported by DOE: DE-SC0019789

We are seeking collaborators who will permit us to measure the plasma density on their experiments